

Please un-install any previous versions of CodeReview software prior to installing this version (which is CodeReview 4.10 software).

After the un-installation process is complete, it is best to use Windows Explorer to completely delete the following folder/files:

C:\Program Files\Pinpoint Technologies, Inc\ZDR

After removing the previous version of CodeReview, install the new CodeReview. by unzipping the "Code410.zip" file, and run the "Code410.exe" file.

The CodeReview program requires an OS of Windows 2000 Pro (or above), as well as an InfraRed Reader Device.

Please note that CodeReview is not designed to operate with Windows VISTA

FYI, there is a USB-style of InfraRed Reader Device that can be purchased through ZOLL Canada's Customer Support department, if one is required for your site.

Give us a call if additional information is desired.

- John Morris (1 866 442-1011 ext 223)

Technical Support

ZOLL Medical Canada
1750 Sismet Road - Unit 1
Mississauga, Ontario
L4W 1R6

FAX: 905-629-0575

CanadaTechSupport@zoll.com